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Lagerbedingung:

Product description:

ORBIS° Duplicating-silicones are universal, duplicating materials for all kinds of duplications. Due to the practical way of mixing, the superior reproduction of details, the constancy of shape and long shelf life, duplicating with ORBIS° Duplicating-silicone is an extremely rational and time-saving method. Without further treatment the surface of the models produced in this duplication form is smooth and hard. Due to its platinum hardener system, the materials are virtually shrink less and it is not necessary to pour out the duplicating form immediately.

The final hardness of the cured silicone is approx. 20 ± 2 shore A.

Dosage:

ORBIS° Duplicating-silicones are mixed in a ratio of 1:1 of the components A + B.

Processing:

The mixing of the two components can be effected either in a beaker or with the help of the numerous dosing and mixing devices available. The most rational and cleanest method is mixing in a machine. These units doses, mixes and dispenses automatically and without any air bubbles at the touch of a button.

Due to the short processing time of the "speed version" a manual mixing should be limited to a max. quantity of 1 duplication.

	Duplicating silicone Type A and B	Duplicating silicone <i>speed</i> Type A and B
Mixing time	0:45	0:30
Working time (incl. Mixing)	6:00 +/- 1:00	2:30 +/- 0:30
Deflasking	25:00 +/- 5:00	10:00 +/- 2:00
Model production after	30:00	15:00

Mixing- and processing time refer to a room temperature of 73°F (23°C) and a relative air humidity of 50 %.

Preparatory steps:

Almost every material, e. g. plaster models, metals or PMMA can be duplicated with ORBIS° Duplicating-silicones. The master model should not be soaked with water but should be dry and clean of all insulating agents. To save material, the base part of the model can be embedded in a foamed plastic ring supplied with the glass- and quick flask in order to cover only the teeth and jaw parts to be duplicated.

Duplication:

If the material is mixed manually, pour mass into the flask in a thin stream. Using a machine the stream is determined by this unit. The distance should be approx. 20 cm

Further processing:

The best results regarding precision and homogeneity are achieved by means of the so called pressure duplicating technique. If the duplicating form was produced under pressure it is decisive that the model is produced under the same pressure, too. The master model may be deflasked directly after termination of the vulcanisation. The form is maintained and can be poured out several times. Before pouring out, the duplicating form is to be sprayed with surface pressure equalizing agent. By this the tension of the surface is relieved and air bubbles which could occur during pouring out can be reduced. All kinds of dental plaster, model materials and embedding masses are suitable for this method.

General information:

ORBIS° Duplicating-silicones shows no negative effects in the physiological and toxicological respect. Chemical and thermal disinfection is possible.

Storage conditions: